
MIThora™
Instructions for use
ICD050 MIThora IFU Rev 1.0

CONCEPT

DESIGN

PRACTICE

Caution: This device is restricted to sale by or on the order of a physician.

1. DESCRIPTION


The MIThora is a reusable trocar designed to establish a port of entry for thoracoscopic instruments used during minimally invasive surgery. It is delivered nonsterile, to be sterilized and reprocessed by moist sterilization (autoclaving) by the user.

The MIThora has a trocar cannula and a one trocar obturator equipped with a ball handle. The trocar cannula has an internal diameter of 17mm and is 127mm long.

The MIThora can be delivered in a polybag, or within a sterilization case.

Table 1: Components of MIThora MIT0000000

| Name of components | Material | Part No. |
|-----------------------------------|--|-----------------|
| MITHora Trocar Obturator | Stainless Steel 17-4PH ASTM A564 and Silicone ASTM D2000 | MIT0000002 |
| MITHora Trocar Cannula | Stainless Steel 17-4PH ASTM A564 & ASTM A693 | MIT0000005 |
| Name of Sterilization Case | Material | Part No. |
| MITHora Case | 5052 H32 ASTM B209 aluminium Anodized MIL-A-8625F Type II Class III, Stainless steel 304 2B ASTM A666, Nylon ASTM D3359, Silicone ASTM D2000, Polypropylene ASTM D4101, R-5000 Radel ASTM D6394 SP0312 | MIT0000006 |

2. MATERIAL

The ball handle is the only part of the device containing silicone ASTM D2000 and stainless-steel ASTM A564. All other parts are made of stainless-steel ASTM A564 or ASTM A693 only.

3. INTENDED USE

The MITHora is intended to establish a port of entry for thoracoscopic instruments used during minimally invasive spinal surgery using an anterior approach. The MITHora is intended to be inserted through the rib cage between the ribs, after one lung has been deflated and a small incision has been performed between two consecutive ribs on the side of the deflated lung.

The MITHora is **not** intended to maintain internal rib cage positive pressure during surgery with the use of CO₂ insufflation and should only be used when the maintenance of pressure in the rib cage is already compromised or unnecessary.

The MITHora is an external communicating device in contact with tissues and bone for limited contact duration (less than 24h).

4. CONTRAINDICATIONS FOR USE

Endoscopic trocars are contraindicated for any uses generally contraindicated for minimally invasive techniques.

5. ELIGIBLE PATIENT POPULATION

The MIThora is indicated for patients requiring minimally invasive spinal surgery using an anterior approach.

6. WARNINGS AND PRECAUTIONS

WARNING: The MIThora is only intended for use by individuals with adequate training and familiarity with minimally invasive surgery using an anterior approach.

The MIThora™ is only intended for use by individuals familiar with minimally invasive techniques. For further information about minimally invasive techniques, complications and hazards, consult the medical literature.

Prior to use, read and follow these instructions for use as well as those of the instruments to be used during the procedure. Damage to the instrument can lead to patient injuries. Always inspect instrument carefully for overall integrity before use.

Improper use of this product can result in life-threatening injury to internal organs and vasculature. Use extreme caution during trocar insertion.

Always keep the trocar obturator straight relative to the tube when inserting or removing. Introducing or removing the trocar obturator at an angle relative to the cannula can damage the cannula.

Although the trocar obturator is blunt, care must be taken when introducing to avoid damage to major vessels and other anatomic structures.

Keep organs out of reach of trocar penetration by ensuring a suitable positioning of the patient's body.

Direct the trocar away from major vessels and other anatomic structures.

Do not use excessive force.

Instruments with a diameter of 17mm or less can be supported by the MIThora™. The usage of instruments with a diameter higher than 17mm is not possible.

MIThora is radiopaque.

After removing the instruments from the cavity, inspect the surgical site for hemostasis and take appropriate steps to achieve hemostasis as needed.

For incisions made, consider the suture of the underlying fascia at the end of the procedure to reduce the risk for incisional herniation.

7. POTENTIAL ADVERSE EVENTS

| | |
|----------------------------|-------------|
| Superficial lesions | Hematoma |
| Injury to internal vessels | Peritonitis |
| Bleeding | Infection |

8. SAFETY PRECAUTIONS

If using instruments from different manufacturer, verify compatibility of instruments before use to avoid complications during surgery.

Become familiar with the MIThora prior to employing it in a surgical procedure to avoid damage to patient, to operator or to instrument.

Careful handling of instruments is necessary to avoid damage or breakage.

Care should be taken when removing instruments not to prematurely dislodge the cannula.

All precautions applicable to minimally invasive procedures should be observed at all times.

Dispose of the device that has become unusable in accordance with local regulations. Container for sharp items shall be used in case of sharp edges on the damaged device.

9. REPROCESSING & STERILIZATION

| | |
|----------------------------------|---|
| WARNINGS | <p>The MITHora is provided cleaned but non-sterile and must be sterilized prior to use. ANSI/AAMI ST79 recommended practices for in-hospital sterilization should be followed.</p> <p>Reprocessing in healthcare facilities is required to take them from their state after clinical use to the state of being cleaned, sterilized and ready for their next use.</p> <p>Long narrow cannulations, recesses and blind ends require particular attention during cleaning since both chemicals and rinse water may be entrapped in them.</p> <p>Avoid impact, scratching, bending or surface contact with any materials that might affect the instrument surface.</p> |
| Limitations on processing | <p>Repeated processing has minimal effect on this instrument. End of life is normally determined by wear and damage due to use.</p> |

| RECOMMENDATIONS | |
|--|--|
| Pre-Treatment at the point of use | <p>Remove gross soil with disposable low-linting cloth/paper wipe</p> |
| Containment and transportation: | <p>It is recommended that this instrument is processed as soon as is reasonably practical following use.</p> <p>Keep the Sterilization Case lid side up during transport.</p> |
| Preparation before cleaning | <p>Visually inspect for damage or the presence of blood or tissue. If blood or tissue is observed, it must be thoroughly cleaned with an automated washer or manually.</p> <p>Ensure disassembly of the components of the instrument before cleaning.</p> |
| Automated Cleaning | <p>Use an automated washer preferably compliant to ISO 15883-2 with lid on.</p> <p>Load the instruments such that the cannulation of the tube can drain.</p> <p>When unloading check the cannulation of the tube for complete removal of visible soil. If necessary, repeat cycle or use manual cleaning.</p> |
| Manual Cleaning | <p>If automated cleaning unavailable:</p> <p>Rinse gross soil.</p> <p>Using soft brush, apply neutral pH detergent solution to all surfaces.</p> <ul style="list-style-type: none"> ○ Enzyme Soak: Refer to detergent manufacturer instructions for concentration. <p>Rinse using DI water. Ensure that water passes through the cannulation of the tube.</p> |
| Post-cleaning drying | <p>Prior to drying, ensure all liquids are voided by tilting the tube.</p> |

| RECOMMENDATIONS | |
|----------------------------------|--|
| Inspection and Testing | <p>Visually inspect for damage and wear. Edges should be free of nicks and present a continuous edge.</p> <p>Check instruments with long slender features (particularly rotating instruments) for distortion. Check assembly with mating components.</p> |
| Packaging | <p>The MITHora components should be loaded into a dedicated tray.</p> <p>Wrap the tray in accordance with local procedures using standard wrapping techniques such as those described in ANSI/AAMI ST79. Recommended wrapping: Heavy duty CSR paper wrap or Non-Spunbond, Meltblown, Spunbond (SMS) Material.</p> <p>When a tray is not available, the MITHora can be packaged in a sterilization pouch compatible with the steam sterilization process.</p> |
| Sterilization | <p>In a properly functioning calibrated steam sterilizer effective sterilization may be achieved using the following parameters:</p> <p>Cycle: Pre-Vacuum Temperature: 270-275°F (132-135°C) Exposure time: 6 min. Drying time: 60 minutes minimum</p> |
| Post-sterilization drying | <p>Drying within the sterilization vessel is standard practice in hospitals. ANSI/AAMI ST79 provides guidance to hospitals for selecting appropriate drying parameters based on the sterilization cycle that is being conducted.</p> |
| Storage | <p>Store in a controlled environment storage cabinet.</p> |
| Additional Information | <p>When sterilising reusable instruments in one autoclave cycle ensure that the sterilizer's maximum load is not exceeded.</p> |

It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

10. USER INSTRUCTIONS

1. Before beginning the procedure, verify overall compatibility of all instruments and accessories.
2. Inspect the instruments for any damage. Do not use the instrument if any damage is noted. Return the instrument to Spino Modulation Inc. if it is not in acceptable condition for surgery.
3. Select and follow a suitable thoracoscopic protocol.
4. Make a small incision where the instrument will be introduced. A large, deep incision may be necessary for the blunt design of the trocar obturator.

5. Create a secondary incision of adequate size to accommodate the MITHora tube. Note: greater trocar insertion force will be required if the incision is too small. Furthermore, too large of an incision, may increase possible port instability.
6. Insert the trocar obturator and tube assembly through the incision by applying continuous downward pressure until the body cavity has been completely penetrated.
7. Position the cannula as desired.
8. Remove the obturator and insert appropriately sized instruments.
9. Remove the tube and suture the incision site.

11. LIMITED WARRANTY AND DISCLAIMER

SPINO MODULATION PRODUCTS ARE SOLD WITH A LIMITED WARRANTY TO THE ORIGINAL PURCHASER AGAINST DEFECTS IN WORKMANSHIP AND MATERIALS. ANY OTHER EXPRESS OR IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS, ARE HEREBY DISCLAIMED.

IF MORE THAN TWO YEARS HAVE ELAPSED BETWEEN THE DATE OF ISSUE/REVISION OF THIS IFU AND THE DATE OF CONSULTATION, CONTACT SPINO MODULATION FOR CURRENT INFORMATION AT 1 888 988 2747.

12. LABELING SYMBOL


| | |
|---|--|
|  | Lot Code |
|  | Catalogue Number |
|  | Material |
|  | Date of manufacture |
|  | Manufacturer |
|  | Non-sterile |
|  | Caution, consult the user instructions |
|  | Consult instructions for use |

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