
spinomodulation

a subsidiary of spinologics

MIThora™

Instructions for use

ICD050 MIThora IFU Rev 1.1

CONCEPT

DESIGN

PRACTICE

Caution: This device is restricted to sale by or on the order of a physician.

1. DESCRIPTION

The MITHora™ is a reusable trocar designed to establish a port of entry for thoracoscopic instruments used during minimally invasive surgery. It is delivered nonsterile, to be sterilized and reprocessed by moist sterilization (autoclaving) by the user.

The MITHora™ has a trocar cannula and a one trocar obturator equipped with a ball handle. The trocar cannula has an internal diameter of 17mm and is 127mm long.


The MITHora™ is delivered non-sterile in a polybag. The MITHora™ sterilization case is sold separately.

Table 1: Components of MITHora™ MIT0000000

Name of components	Material	Component No.
MITHora™ Trocar Obturator	Stainless Steel as per ASTM A564 and Silicone ASTM D2000	MIT0000002
MITHora™ Trocar Cannula	Stainless Steel as per ASTM A564 & ASTM A693	MIT0000005

Table 2: MITHora™ sterilization case MIT0000006

Name of Sterilization Case	Material	Part No.
MITHora™ Case	5052 H32 ASTM B209 aluminium Anodized MIL-A-8625F Type II Class III, Stainless steel 304 2B ASTM A666, Nylon ASTM D3359, Silicone ASTM D2000, Polypropylene ASTM D4101, R-5000 Radel ASTM D6394 SP0312	MIT0000006

2. MATERIAL

The ball handle is the only part of the device containing silicone ASTM D2000 and stainless-steel ASTM A564. All other parts are made of stainless-steel ASTM A564 or ASTM A693 only.

3. INTENDED USE

The MITHora™ is intended to establish a port of entry for thoracoscopic instruments used during minimally invasive spinal surgery using an anterior approach. The MITHora is intended to be inserted through the rib cage between the ribs, after one lung has been deflated and a small incision has been performed between two consecutive ribs on the side of the deflated lung.

The MITHora™ is **not** intended to maintain internal rib cage positive pressure during surgery with the use of CO₂ insufflation and should only be used when the maintenance of pressure in the rib cage is already compromised or unnecessary.

The MITHora™ is an external communicating device in contact with tissues and bone for limited contact duration (less than 24h).

4. CONTRAINDICATIONS FOR USE

Endoscopic trocars are contraindicated for any uses generally contraindicated for minimally invasive techniques.

5. ELIGIBLE PATIENT POPULATION

The MITHora™ is indicated for patients requiring minimally invasive spinal surgery using an anterior approach.

6. WARNINGS AND PRECAUTIONS

WARNING: The MITHora™ is only intended for use by individuals with adequate training and familiarity with minimally invasive surgery using an anterior approach.

- The MITHora™ is only intended for use by individuals familiar with minimally invasive techniques. For further information about minimally invasive techniques, complications and hazards, consult the medical literature.
- Prior to use, read and follow these instructions for use as well as those of the instruments to be used during the procedure. Damage to the instrument can lead to patient injuries. Always inspect instrument carefully for overall integrity before use.
- Improper use of this product can result in life-threatening injury to internal organs and vasculature. Use extreme caution during trocar insertion.
- Always keep the trocar obturator straight relative to the tube when inserting or removing. Introducing or removing the trocar obturator at an angle relative to the cannula can damage the cannula.
- Although the trocar obturator is blunt, care must be taken when introducing to avoid damage to major vessels and other anatomic structures.
- Keep organs out of reach of trocar penetration by ensuring a suitable positioning of the patient's body.
- Direct the trocar away from major vessels and other anatomic structures.
- Do not use excessive force.
- Instruments with a diameter of 17mm or less can be supported by the MITHora™. The usage of instruments with a diameter higher than 17mm is not possible.
- MITHora™ is radiopaque.
- After removing the instruments from the cavity, inspect the surgical site for hemostasis and take appropriate steps to achieve hemostasis as needed.
- For incisions made, consider the suture of the underlying fascia at the end of the procedure to reduce the risk for incisional herniation.

7. POTENTIAL ADVERSE EVENTS

- Superficial lesions
- Injury to internal vessels
- Bleeding
- Hematoma
- Peritonitis
- Infection

8. SAFETY PRECAUTIONS

- If using instruments from different manufacturer, verify compatibility of instruments before use to avoid complications during surgery.
- Become familiar with the MITHora™ prior to employing it in a surgical procedure to avoid damage to patient, to operator or to instrument.
- Careful handling of instruments is necessary to avoid damage or breakage.
- Care should be taken when removing instruments not to prematurely dislodge the cannula.
- All precautions applicable to minimally invasive procedures should be observed at all times.
- Dispose of the device that has become unusable in accordance with local regulations. Container for sharp items shall be used in case of sharp edges on the damaged device.

9. REPROCESSING & STERILIZATION

WARNINGS	<p>The MITHora™ is provided non-sterile and must be cleaned and sterilized prior to first use and between each surgery. ANSI/AAMI ST79 recommended practices for in-hospital sterilization should be followed.</p> <p>Reprocessing in healthcare facilities is required to take them from their state after clinical use to the state of being cleaned, sterilized and ready for their next use.</p> <p>Long narrow cannulations, recesses and blind ends require particular attention during cleaning since both chemicals and rinse water may be entrapped in them.</p> <p>Avoid impact, scratching, bending or surface contact with any materials that might affect the instrument surface.</p>
Limitations on processing	<p>Repeated processing that includes ultrasonic, mechanical washing and sterilization have minimal effect on MITHora™ system. End of life is normally determined by wear and damage due to use. Evidence of damage and wear on a device may include but is not limited to corrosion (i.e. rust, pitting), discoloration, excessive scratches, flaking, wear and cracks. Improperly functioning devices, devices with unrecognizable markings, missing and removed (buffed off) part numbers, damaged and excessively worn devices should not be used.</p>

RECOMMENDATIONS	
Pre-Treatment at the point of use	Remove gross soil with disposable low-linting cloth/paper wipe
Containment and transportation:	<p>It is recommended that this instrument is processed as soon as is reasonably practical following use.</p> <p>Keep the Sterilization Case lid side up during transport.</p>
Preparation before cleaning	<ul style="list-style-type: none"> • Visually inspect for damage or the presence of blood or tissue. If blood or tissue is observed, it must be thoroughly cleaned with an automated washer or manually. • Ensure disassembly of the components of the instrument before cleaning.
Automated Cleaning Process	<ul style="list-style-type: none"> • At point of use take appropriate measures to ensure contaminant does not dry on soiled device. • Ensure disassembly of the components of the instrument before cleaning. • Pre-rinse the system contents with cool water for a minimum of 3 minutes. • Use an automated washer compliant to ISO 15883-2 with lid on to contain implant and instruments components.

RECOMMENDATIONS																																	
	<ul style="list-style-type: none"> • Load the instruments such that cannulations can drain. • When unloading check cannulations for complete removal of visible soil. If necessary, repeat cycle or use semi-automated cleaning. <p>Automated Cleaning with an Enzymatic Solution:</p> <ol style="list-style-type: none"> 1. Execute the cycle using a pH neutral enzymatic detergent according to the following parameters: <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th colspan="4">Motor Speed High</th> </tr> <tr> <th>Step</th> <th>Description</th> <th>Minimum Temperature</th> <th>Minimum Cycle</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>Pre-Wash</td> <td>Cold Tap Water (Facility)</td> <td>2 Minutes</td> </tr> <tr> <td>2</td> <td>Enzyme Wash</td> <td>Hot Tap Water (Facility)</td> <td>4 minutes</td> </tr> <tr> <td>3</td> <td>Wash (Detergent % according to manufacturer specification)</td> <td>65.5 °C (149.9°F)</td> <td>2 Minutes</td> </tr> <tr> <td>4</td> <td>Rinse</td> <td>Hot Tap Water</td> <td>10 Minutes</td> </tr> <tr> <td>5</td> <td>Thermal Rinse with Lubricant (% Lubricant according to manufacturer specification)</td> <td>82.2 °C (180.0°F)</td> <td>10 Minutes</td> </tr> <tr> <td>6</td> <td>Hot Air Dry</td> <td>HIGH</td> <td>6 Minutes</td> </tr> </tbody> </table> <ol style="list-style-type: none"> 2. Upon completion, unload the washer-disinfector. 3. Visually inspect the Devices for remaining soil. 	Motor Speed High				Step	Description	Minimum Temperature	Minimum Cycle	1	Pre-Wash	Cold Tap Water (Facility)	2 Minutes	2	Enzyme Wash	Hot Tap Water (Facility)	4 minutes	3	Wash (Detergent % according to manufacturer specification)	65.5 °C (149.9°F)	2 Minutes	4	Rinse	Hot Tap Water	10 Minutes	5	Thermal Rinse with Lubricant (% Lubricant according to manufacturer specification)	82.2 °C (180.0°F)	10 Minutes	6	Hot Air Dry	HIGH	6 Minutes
Motor Speed High																																	
Step	Description	Minimum Temperature	Minimum Cycle																														
1	Pre-Wash	Cold Tap Water (Facility)	2 Minutes																														
2	Enzyme Wash	Hot Tap Water (Facility)	4 minutes																														
3	Wash (Detergent % according to manufacturer specification)	65.5 °C (149.9°F)	2 Minutes																														
4	Rinse	Hot Tap Water	10 Minutes																														
5	Thermal Rinse with Lubricant (% Lubricant according to manufacturer specification)	82.2 °C (180.0°F)	10 Minutes																														
6	Hot Air Dry	HIGH	6 Minutes																														
Semi-Automated Cleaning Process	<p>If automated cleaning unavailable:</p> <ul style="list-style-type: none"> • At point of use take appropriate measures to ensure contaminant does not dry on soiled device. • Ensure disassembly of the components of the instrument before cleaning. 																																

RECOMMENDATIONS	
	<ul style="list-style-type: none"> • Rinse gross soil from implants and instrument using warm (25-40°C) utility water. Rinse each component of the system for a minimum of 1 minute. • Enzyme Soak: Refer to detergent manufacturer instructions for concentration and temperature. Soak for a minimum of 1-5 minutes. • Using a soft brush, scrub all surfaces, ensure instruments are cleaned. • Sonication: Completely submerge devices in the sonicator solution. Space the devices so there is no contact between the instruments; do not overload the ultrasonic cleaning unit. Use the processing time recommended by the manufacturer of the ultrasonic unit. • Rinse using DI water for a minimum of 10 minutes. Ensure water passes through cannulations. • Prior to drying, ensure all liquids are voided by tilting the instruments. • Dry each device with a low-lint wipe or cloth, or with compressed air until visually dry. • All instruments: Visually inspect for visible soil, damage, and wear. If visible soil is present, repeat cleaning process. Cutting edges should be free of nicks and present a continuous edge. <p>Check instruments with long slender features for distortion. Check assembly with mating components. Apply lubricant in accordance with lubricant manufacturer instructions.</p>
Sterilization Process	<ul style="list-style-type: none"> • Devices should be completely dry prior to wrapping and sterilization. • Instruments should be loaded into dedicated tray. • Wrap the tray in accordance with local procedures using standard wrapping techniques such as those described in ANSI/AAMI ST79. • Wrap with heavy duty sterilization wrap. (Recommended wrapping: Heavy duty CSR paper wrap or Non-Spunbond, Meltblown, Spunbond (SMS) Material.) • Process devices in a properly functioning calibrated steam sterilizer according to the following parameters: <p>Prevacuum (4-pulse) conditioning; 4-minute, 132 °C exposure; 30-minute dry time</p>

RECOMMENDATIONS	
Storage	Store in a controlled environment storage cabinet.
Additional Information	When sterilizing reusable instruments in one autoclave cycle ensure that the sterilizer's maximum load is not exceeded.

It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

10. USER INSTRUCTIONS

1. Before beginning the procedure, verify overall compatibility of all instruments and accessories.
2. Inspect the instruments for any damage. Do not use the instrument if any damage is noted. Return the instrument to Spino Modulation Inc. if it is not in acceptable condition for surgery.
3. Select and follow a suitable thoracoscopic protocol.
4. Make a small incision where the instrument will be introduced. A large, deep incision may be necessary for the blunt design of the trocar obturator.
5. Create a secondary incision of adequate size to accommodate the MITHora™ tube. Note: greater trocar insertion force will be required if the incision is too small. Furthermore, too large of an incision, may increase possible port instability.
6. Insert the trocar obturator and tube assembly through the incision by applying continuous downward pressure until the body cavity has been completely penetrated.
7. Position the cannula as desired.
8. Remove the obturator and insert appropriately sized instruments.
9. Remove the tube and suture the incision site.

11. LIMITED WARRANTY AND DISCLAIMER

SPINO MODULATION PRODUCTS ARE SOLD WITH A LIMITED WARRANTY TO THE ORIGINAL PURCHASER AGAINST DEFECTS IN WORKMANSHIP AND MATERIALS. ANY OTHER EXPRESS OR IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS, ARE HEREBY DISCLAIMED.

IF MORE THAN TWO YEARS HAVE ELAPSED BETWEEN THE DATE OF ISSUE/REVISION OF THIS IFU AND THE DATE OF CONSULTATION, CONTACT SPINO MODULATION FOR CURRENT INFORMATION AT 1 888 988 2747.

12. LABELING SYMBOL

	Lot Code
	Catalogue Number
	Material
	Date of manufacture
	Manufacturer
	Non-sterile
	Caution, consult the user instructions
	Consult instructions for use

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Spino Modulation, Inc.